



Republic of the Philippines
DEPARTMENT OF HEALTH
PHILIPPINE CHILDREN'S MEDICAL CENTER

Quezon Avenue, Quezon City 1100

website: www.pcmc.gov.ph email: officeofthedirector@pcmc.gov.ph

Trunkline: 8588-9900 DirectLine: 8924-0836 Fax No: 8924-0840

BID BULLETIN NO. BB-2022-015

Date : **March 14, 2022**

TO : **ALL PROSPECTIVE BIDDERS**

FROM : **BIDS AND AWARDS COMMITTEE**

SUBJECT : **AMENDMENT AS STATED**

Name of Project : **Supply and Delivery of Forty (40) unit of Infusion Pump
(ABC = Php 2,680,000.00) IB 2022-059**

This Bid Bulletin No. **BB-2022-015** outlines the amendments in the bidding documents for above projects as follows:

1. Under Section VII. Technical Specifications

FROM	TO
A. Specifications:	A. Specifications:
1. Can be used for any kind of IV set (universally compatible no need for recalibration), and in any solution-crystalloid, colloid TPN	1. Can be used for any kind of IV set (universally compatible), and in any solution-crystalloid, colloid TPN
3. Capable of delivering Manual and Automatic three (3) boluses including Target Volume and Target Time	3. Capable of delivering Manual and Automatic one(1) to three (3) boluses including Target Volume and Target Time
4. Keypad lock or parameter lock, locks the keypad and prevents unauthorized access to safety relevant functions	4. Keypad lock/ Screen lock or parameter lock, locks the keypad and prevents unauthorized access to safety relevant functions
5. Deactivation drop sensor eliminates unnecessary drop sensor alarms during rugged transport of patients	5. Deactivation drop/ drip sensor eliminates unnecessary drop sensor alarms during rugged transport of patients
6. Viewable battery charge level in terms of hours and minutes, displays how long the pump can operate in battery mode	6. Viewable battery charge level/ status on screen in terms of hours and minutes, displays how long the pump can operate in battery mode or Manufacturers classification
19. Alarm in case of incorrect dosage and apparatus malfunction the pump switches off automatically	19. Alarm prompts as high level alarm in case of incorrect dosage and apparatus malfunction the pump stops off automatically

For information and guidance of all concerned.

MARIA ROSARIO S. CRUZ, MD

Chairperson, Bids and Awards Committee

