

Republic of the Philippines DEPARTMENT OF HEALTH PHILIPPINE CHILDREN'S MEDICAL CENTER

Quezon Avenue, Quezon City 1100 website: <u>www.pcmc.gov.ph</u> email: <u>officeofthedirector@pcmc.gov.ph</u> Trunk Line: 8588-9900 to 20 Direct Line: 8924-6601

NOTICE OF AWARD

NOA-2024-176

A.M.I EQUIPMENT SERVICES AND SOLUTIONS INC.

U-1012 Medical Plaza Bldg. #25 San Miguel Avenue, Ortigas Center Brgy.San Antonio Pasig City E-mail Add.: mattabay@gmail.com

Dear Sir / Madam:

This is to inform you that based on the result of the conducted Competitive Bidding for the project stated below under Invitation to Bid No. IB-2024-091 as per BAC Resolution No. R2024-12-708, your proposal was found to be responsive.

ITEM DESCRIPTION	TOTAL COST
Supply, Delivery, Installation, Testing and Commissioning of Washer Disinfector (Turnkey) - (see attached Annex A for the detailed offer)	Php7,999,000.00

You are hereby required to provide on or before either of the following form:

D 6 JAN 2025

the Performance Security in

FORM OF PERFORMANCE SECURITY	AMOUNT OF PERFORMANCE SECURITY (Equal to Percentage of the Total Contract Price)
 a) Cash or cashier's/ manager's check issued by a Universal or Commercial Bank b.) Bank draft/guarantee or irrevocable letter of credit issued by a Universal or Commercial Bank: Provided, however, that it shall be confirmed or authenticated by a Universal or Commercial Bank, if issued by a foreign bank 	Five percent (5%) <u>Php399,950.00</u>
c) Surety bond callable upon demand issued by a surety or insurance company duly certified by the Insurance Commission as authorized to issue such security.	Thirty percent (30%) Php2,399,700.00

Kindly disregard and report any form of solicitation using the name of the Executive Director and/or the PCMC. These unsanctioned requests are unlawful and will not be tolerated





NOTICE OF AWARD NOA-2024-176 A.M.I EQUIPMENT SERVICES AND SOLUTIONS INC.

Failure to provide the Performance Security shall constitute sufficient ground for cancellation of the award and imposition of penalties/ forfeiture corresponding to the Bid Security posted.

Very truly yours,

SONIA B. GONZALEZ, MD, MScHSM, MPM Executive Director

Conforme:

This is to certify that the company has authorized me to accept this award, sign all related documents and hold the company bound by rules and laws applicable thereto.

Authorized Signatory (Signature over printed name)

Designation

Date

Kindly disregard and report any form of solicitation using the name of the Executive Director and/or the PCMC. These unsanctioned requests are unlawful and will not be tolerated

-	UNIT		_	ITEM DESCRIPTION
1	Unit			Delivery, Installation, Testing and Commissioning of Washer Disinfector (TurnKey <i>d: MMM</i>
				a: MMM 2/Model: Uniclean PL II 10
				Amount Awarded = Php 7,999,000.00
		-	The second s	cifications of the Main Equipment
			1.	Chamber material: stainless steel AISI 316L with chamber capacity of 216 liters ,
			2.	Exterior Dimensions offered: 700mm Width x 2,300mm height x 777mm depth
			3.	Door configuration:
				Automatic vertical sliding door
				Double panel tempered glass
			4	Loading height: adjustable 755mm to 1000mm
			5	Stainless steel pump:
				• 3hp ,
				• 2-3 chemical dosing pump
				• Rotary spray arms
				•No compressed air required
			6	Control Screen panel: MMM Smart HMI with 7" color screen at eye level with software system for remote maintenance
			7	Rack loading access from all sides
			8	Chamber lighting 3 LED colored inside indicators
				• Blue - Standby
				• White - Ongoing Cycle
				• Red - Errors
			9	Cycle time minimum of 30mins including drying time
			10	3 to 5 level racks
			11	Build-in thermal printer that can validate cycles procedures .
			12	Can be connected to existing reverse osmosis
			13	HEPA 13 filter for air Discharge
_			14	Provision of separate / CKT power supply line intended for the unit
				Operating at 220-420 volts 3-phase,50–60 Hz
				Power consumption: 6-9 kw
				All mechanical parts are heavy duty and made up of stainless steel 3161
				Maintainable and serviceable
		в. Ас	cess	ories to be provided per unit
			1. 5	Stainless steel 15 trays
			2.	Cart - 2pcs stainless steel trolley
				3inch four caster wheel
				Easy glide trasport rollers for loading Loading rack accessible from all side

٠

QTY	UNIT			ITEM DESCRIPTION
1	Unit	Sup	ply,	Delivery, Installation, Testing and Commissioning of Washer Disinfector (TurnKey
				Rack can be placed under the trasport cart \sim
			3.	Safety Electrical Devices with the rated capacity of the unit
				> AVR
				Brand: Panther
				Model: PVE - 1000
		C. C	Consi	umables to ber Provided per unit
			1	7 gallon enzymatic cleanser
			2	2 liters rinse and drying activator
			3	100 pcs. Test indicator
			4	10 Sets/pairs of thick rubberized gloves.
			5	4pcs thermal paper
		D	-	
		D.	-	ranty
				With three (3) years warranty on parts and services
			And the second sec	Quarterly Preventive Maintenance Services and Calibration during the warranty period.
		E. 5	and the second s	e of Works
				Conformed to the attached Scope of Works
		F. W	procession and an owner of	luly notarized certificate from Bidder - A.M.I Equipment Services & Solutions
				ifies that:
			i.	Parts, accessories, and consumables are readily available at the authorized Philippine service
				center for a period of five (5) years after the warranty period;
			ii.	It has available competent in-house technical specialists in handling and providing technical
				support as well as maintenance of the equipment being offered;
			iii.	It will conduct training for proper operation and maintenance to the End-users of the equipment
				upon delivery; and
			and the second sec	It will provide replacement/back-up unit while the delivered unit is being repaired.
				luly notarized certificate from Principal Manufacturer - MMM Munchener Medizin
		Mec	hanio	c GmbH
			Statement in call of the local division of the	ifies that:
			i.	All the terms and conditions stated in the bidding documents per IRR of RA 9184 and the
				corresponding contract for the contract for the project shall be honored by the Principal
				Manufacturer, including in the event that a change of dealership will occur during the duration of
				the contract up to the warranty and preventive maintenance period (2 years)
				The expected useful life of the equipment is 10-15 years under 24hours continuos usage in a
			ii.	standard hospital patient room of approximate 50sqr. Meters
				standard nospital patient room of approximate 50sql. Meters
			iii.	Guarantee on the availability of all spare parts, accessories and consumables at least for the next
				five (5) years from testing, commissioning, acceptance and delivery;
			iv.	It has competence in handling and providing technical support as well as maintenance of the
				equipment being offered; and
			v.	Proper disposal of Washer Disinfector
			*•	
				1. After the hospital has decided to disinvest on the Washer Disinfector it is a condemned unit
				2. Decontamination will be done to the decommissioned unit
				3. All connections (electrical, pipes and fittings) will be disconnected and the area where the
				condemned unit was once situated will be prepared for future installation of either the same item
				or other instrument reprocessing item
				4. It will depend on the hospital if it will be recycled or refurbished to be donated or sold to other
				institutions.
				5. If not, the hospital will give the authorizquation to dismantle the unit per part and then dispose of
1				
				lit by colling it of donon like the your thing that other institutions do the conchility for constitution
				it by selling it as scrap like the usual thing that other institutions do the capability for corrective and preventive maintenance of the unit within the warranty period.

A

QTY	UNIT		ITEM DESCRIPTION	
1	Unit	Sup	oly, I	Delivery, Installation, Testing and Commissioning of Washer Disinfector (TurnKey)
		H. D	elive	ry Period:
			1.	Within Sixty to Ninety (60 to 90) Calendar days upon receipt of Purchase Order/Notice to Proceed

Conforme:

5 5 . . .

Authorized Signatory (Signature over printed name)

Designation/Position

Date_

PHILIPPINE CHILDRENS MEDICAL CENTER Quezon Avenue, Quezon City

Name of Office: NURSING SERVICE-STERILIZATION ROOM Description of the Project: SUPPLY, DELIVERY, INSTALLATION, AND COMMISSIONING OF WASHER DISINFECTOR (TURNKEY)

ABC: **8,000,000.00** Reference: **PR No. SR-2024-06** Date: **October 2, 2024**

, ,

SCOPE OF WORKS

INSTALLATION OF WASHER DISINFECTOR

- Positioning and Mounting:
- Position the washer disinfector as per design drawings, ensuring the unit is level and aligned for proper operation.
- · Secure the washer disinfector to the floor if required.

UTILITY CONNECTIONS:

- Water Supply: Connect the washer disinfector to the designated water supply lines, ensuring proper backflow prevention systems.
- Drainage: Connect the drainage outlet to the facility's drainage system, ensuring compliance with local plumbing codes.
- Electrical, connecting from the machine to the toping point for a dedicated power source as per manufacturer specifications

ACCESSORIES INSTALLATION

Accessories for Decontamination

- · Hanging shelves for the decontamination area
- Stainless steel type 304-gauge 20 hairline finish
- · Provision of wall shelf H80cm x L249cm x W55cm divided into two equal horizontal shelves
- · Provision and installation of a sliding window located in the decontamination area
- · double tempered glass horizontal slider window
- size dimension L74cm x W10cm x H50cm
- · Provision of the base cart of ultrasonic cleaner 20L
- stainless steel type 304-gauge 20 hairline finish
- size dimension: L61cm x W38cm x H92cm
- · approximately 2-3 inches four rigid caster wheels
- Install any additional components such as special trays etc. according to the requirements of the unit and the instruments being cleaned.

TESTING AND COMMISSIONING

• Verify that all utility connections are functioning correctly (water flow, drainage, electricity, etc.).

A.M.I. EQUIPMENT SERVICES AND SOLUTIONS INC. 3

- Perform initial run cycles to ensure proper operation, including correct washing, disinfection, and drying functions.
- Ensure the washer disinfector is calibrated according to manufacturer guidelines to meet operational standards for temperature, cycle times, and detergent usage.
- Verify safety features such as door interlocks, emergency stops, and pressure relief systems are functioning properly.

OPERATIONAL TRAINING

•

- Provide comprehensive training to the staff on how to properly operate the washer-disinfector, including loading, program selection, and safety protocols.
- Train staff on routine maintenance tasks such as cleaning filters, replacing detergent containers, and basic troubleshooting.

CONFORME WINFTYED DIZAN

Authorized Signatory Signature over printed name

AMI EDMPMENT SELVICES AND SOLUTIONS, INC

bidding a amicorp. com. ph

0287067300

Contact No:

0287067-300

Name of Company/Firm -No. Company's Official Email Address. Company's Official Contact

(where notices will be sent).

A.M.I.	EQUIP	LANA JELI	C. San
CERT	FL	n	OPY
E.S.	- h	10	
Close	1	Philippine	same ,